

# 40Hz-Auditory Steady State Response

Auditory Evoked Potential, a Powerful Biomarker for Information Processing

## A Powerful Window into Information Processing

The main challenge in preclinical research is finding biomarkers that can be applied across both **clinical and preclinical** settings. The 40Hz-ASSR is a non-invasive biomarker detectable in both humans and rodents. It measures how well the auditory cortex **processes information** and generates synchronous activity at a specific frequency. Disruptions in ASSR are linked to disorders such as **schizophrenia, bipolar disorder,** and **autism spectrum disorder**.

## Why EEG Phenotyping ?



### Translational

This Evoked Potential can be observed in **humans and animals** at a peak frequency of 40Hz. Several **non-invasive** experimental protocols are commonly **used in the clinic** on human patients.



### Objective

40Hz-ASSR provides a **robust** and **specific** measure of brain resonance to auditory stimuli in **freely-moving** animals.



### Sensitivity to Neuropsychiatric Conditions

The 40-Hz ASSR is **sensitive to alterations** in neural circuits that are often disturbed in **neuropsychiatric conditions**. This sensitivity enhances its utility in both understanding and diagnosing these conditions.




## How Do We Produce an Auditory Response?

### ASSR Paradigm



Animals will be stimulated with a train of sound pulses of 5ms at 40Hz, lasting 2s with a resting period of 8s. This procedure is repeated 360 times, lasting 1 hour in total.


### Resulting Cortical Response



First stimulus triggers a greater response, followed by a train of resonance.

# Case Study: Characterizing a Neurodevelopmental Disorder Model

**Step 1**



Electrode **implantation**  
Model selection (EEG)

**Step 2**



40Hz-ASSR recording


**Step 3**



Analysis : Evoked Power and  
Inter-Trial Coherence

## Results


### Grand Average of Responses




**Figure 1** - Close-up of the grand average of responses to the 40Hz-ASSR for both genotypes. **Dark blue** : WT animals (n=15); **Light blue** : Transgenic animals (n=15).

The amplitude of response to the 40Hz ASSR stimulation is significantly decreased for the transgenic animals.

### Evoked Power & Inter-Trial Coherence



**Figure 2** - Evoked power measured in WT and Transgenic rats for 2 recordings, n=15 animals per genotype.



**Figure 3** - Inter-Trial Coherence index measured in WT and Transgenic animals for 2 recordings, n=15 animals per genotype.